

Easy to Use

High Performance

Digital Display


The Digital Hand-Held "Pocket" Refractometer PAL-RI, can easily measure samples and display readings digitally. And temperature of the measurement will be displayed continuously.

Digital Hand-held "Pocket" Refractometer
Cat.No.3850

PAL-RI

Example of display

i.e.) display of the measurement values "1.3425" and "20°C".



The numeric characters will stream from right to left on the display in intervals of 0.8 seconds.

Quick & Easy cleanup!



Actual size

Features

- Requires only a few drops of sample
- Measurement in 3 seconds
- Ergonomically designed for easy, one-handed operation
- Light & Compact, 100g
- Designed with a strap holder option
- Highly water resistant (IP65)
- Convenient storage case
- Battery power indicator

Specifications

Measurement range	: Refractive index 1.3306~1.5284
Resolution	: Refractive index 0.0001
Measurement accuracy	: Refractive index ±0.0003 (water at 20°C)
Measurement temperature	: 5~45°C
Ambient temperature	: 10~40°C
Power supply	: AAA alkaline battery × 2
Battery life	: Approx. 11,000 measurements
International Protection Class	: IP65 Water resistant
Dimensions and weight	: 55(W)×31(D)×109(H)mm, 100g (Main Unit only)

* The displayed refractive index is not the temperature compensated value. It is the refractive index for the displayed measurement temperature.

All ATAGO refractometers are designed and manufactured in Japan.



ATAGO products are acknowledged by HACCP, G.M.P. and G.L.P.



※Specifications and appearance are subject to change without notice.

labo and Co
www.laboandco.com
N° Indigo 08 20 20 16 16