


SOLUTIONS TECHNIQUES pour

laboratoires
& industries



LABO AND CO

Nos compétences et services associés



Une équipe de spécialistes

COOLSAFE PRO 55-4

SPECIFICATIONS



CoolSafe Pro 55-4	Unit	Part no: 7001100058
Digital read-out of temperature		Yes
Digital read-out of pressure		Yes
Ultimate temperature at 20°C room temperature	°C	-55
Power requirements	V/Hz	230/50-60
Power consumption	W	500
Insulation	mm	90
Condenser capacity	kg	3
Condenser capacity per 24 hours	kg	2.5
Total volume of Condenser	L	4
Weight net	kg	40
Cabinet dimensions (DxWxH)	mm	500x400x520
Cabinet material		Polyester coated steel
Condenser material		Stainless steel AISI 316
Cooling media		R 507
Vacuum Valve		Solenoid
Drain tap		Yes
Keypad for control options		No
R 485/USB port for Labogene software		No

Revision date | rev. | version: May 2020 | 0.2 | CoolSafe Pro 55-4 - 230V



LABO AND CO

50 avenue de Grosbois
94440 Marolles-en-Brie

Tel : 01 45 98 74 80

Fax : 01 45 98 77 23

contact@laboandco.com



LABO AND CO



Gardez un œil sur laboandco.com

laboandco.com