



# Vortex Mixers | Accessories



### VT1.2 Tube holding rod

used with tube adapters

**Cat. No.** 18900044

### VT1.3 Universal top plate

Ø 100mm

**Cat. No.** 18900035



### VT1.3.1 tube adapter

for 48 holes test tubes, Ø6mm

**Cat. No.** 18900020



### VT1.3.2 tube adapter

for 15 holes test tubes, Ø10mm

**Cat. No.** 18900021



### VT1.3.3 tube adapter

for 12 holes test tubes, Ø12mm

**Cat. No.** 18900022



### VT1.3.4 tube adapter

for 8 holes test tubes, Ø16mm

**Cat. No.** 18900023



### VT1.3.5 tube adapter

for 8 holes test tubes, Ø20mm

**Cat. No.** 18900024



### VT1.3.6 Platform pad

for <math>\varnothing</math>99mm tubes and small vessels

**Cat. No.** 18900043



### Application of accessories

Accessories	Adjustable speed model (0~2500rpm)	
	Touch mode (High speed area)	Continuous mode (Low speed area)
VT1.1	Y	Y
VT1.2 + VT1.3.1	Y	
VT1.2 + VT1.3.2	Y	
VT1.2 + VT1.3.3	Y	
VT1.2 + VT1.3.4	Y	
VT1.2 + VT1.3.5	Y	
VT1.3 + VT1.3.6	Y	
VT1.3 + VT1.3.1		Y
VT1.3 + VT1.3.2		Y
VT1.3 + VT1.3.3		Y
VT1.3 + VT1.3.4		Y
VT1.3 + VT1.3.5		Y



### VT1.3.7 vacuum chuck

made of rubber

**Cat. No.** 18900158

