

[Click Below](#)

Bottle-top Dispensers

DispensMate/dFlow



Your Accessible Solution

For General Use

Manual pipettes
Pip. Filler/ Stepper
Electronic Dispenser
Magnetic Stirrers
Shakers

For Analytical Labs

Electronic Titrator
Manual Dispenser
Rotary Evaporator
Hotplate
Overhead Stirrers

For Life Science Labs

Electronic Pipettes
SafeVac
Centrifuges-LS
Mixers / Vortex
Thermo Mixers

For Clinical Labs

Centrifuges-Clinical
Rollers and Rotators
Thermo Cyclers
Spectrophotometers



Easy to control

DispensMate

dFlow

DispensMate

Advantages for the user	Features
Excellent chemical resistance	Components are made of PTFE, FEP, BSG, PP
High applicability	Covers a volume range from 0.5 to 5ml
Fully Autoclavable	Fully Autoclavable at 121°C

dFlow

Advantages for the user	Features
Precise control	Motor operation minimizes the grip force
Easy to operate	Remote control panel with large display without the manual error
Two dispensing modes	Dispenser and stepper function options
Excellent chemical resistance	Components are made of PTFE, FEP, BSG.

Click Below

For General Use

Manual pipettes

Pip. Filler/ Stepper

Electronic Dispenser

Magnetic Stirrers

Shakers

For Analytical Labs

Electronic Titrator

Manual Dispenser

Rotary Evaporator

Hotplate

Overhead Stirrers

For Life Science Labs

Electronic Pipettes

SafeVac

Centrifuges-LS

Mixers / Vortex

Thermo Mixers

For Clinical Labs

Centrifuges-Clinical

Rollers and Rotators

Thermo Cyclers

Spectrophotometers



LABO AND CO
your experience®

Bottle-Top Dispensers

DispensMate/dFlow



Your Accessible Solution

Technical data

DispensMate

Volume Range	Graduation	A \leq \pm		CV \leq	
ml	ml	%	μ l	%	μ l
0.5-5	0.1	0.5	25	0.1	5
1.0-10	0.2	0.5	50	0.1	10
2.5-25	0.5	0.5	125	0.1	25
5.0-50	1	0.5	250	0.1	50

dFlow

Volume range	0.01mL-99.99mL, Max piston lift is10mL, resolution is
Volume accuracy	Dispenser R=0.6% CV=0.2%
	Stepper R=0.5% CV=0.3%
Operating temperature	10 ° ~30 °
Quality control	DIN EN ISO 8655
Control type	External control

Ordering Information

DispensMate

Cat. No.	Volume
7032100001	0.5-5ml
7032100002	1-10ml
7032100003	2.5-25ml
7032100004	5-50ml
17400037	Reagent Bottle (Brown, 1L)

dFlow

Cat. No.	Volume
7017100001	Digital Burette 0.01-99.99ml

Accessories

Cat. No.	Description
17000069	S40, Adapter, 45/40 mm
17400017	GL32, Adapter, 45/32 mm
17400018	GL38, Adapter, 45/38 mm
17400019	GL25, Adapter, 32/25 mm
17400020	GL28, Adapter, 32/28 mm
17400021	Filling tube (m)
17100073	Discharge tube (m)