

[Click Below](#)

For General Use

Manual pipettes
Pip. Filler/ Stepper
Electronic Dispenser
Magnetic Stirrers
Shakers

For Analytical Labs

Electronic Titrator
Manual Dispenser
Rotary Evaporator
Hotplate
Overhead Stirrers

For Life Science Labs

Electronic Pipettes
SafeVac
Centrifuges-LS
Mixers / Vortex
Thermo Mixers

For Clinical Labs

Centrifuges-Clinical
Rollers and Rotators
Thermo Cyclers
Spectrophotometers



Manual Pipette Filler Levo Plus/Levo ME/Levo

DVAB
Your Accessible Solution



Levo Plus

Advantages for the user	Features
Easy to operate	Light weight and single-handed operation
High applicability	For most plastic and glass pipettes from 0.1 to 100ml
Selective speed	Eight speeds of selection for aspiration and dispensing.
Friendly display	Large LCD display showing low battery warning and speed setting.
Rechargeable battery	Powerful rechargeable Li-ion battery enables 8h runtime
0.45µL Filter	Universally used filter and is autoclavable

Technical data

Aspirate Speeds	25ml < 5s (8 shift)
Dispense Speeds	Motor (8 shift) / Gravity
Battery	Lithium-ion
Battery Service Life	More than 8 Hours of Intermittence Use
Charging Time	2~3 Hours
Pipette Types	Glass or Plastic Pipette (0.1-100mL), Pasteur Pipettes
Filter	0.45 µm Hydrophobic
Weight	200g

Levo

Advantages for the user	Features
Easy to operate	Lightweight and single-handed operation
High applicability	For most plastic and glass pipettes from 0.1 to 100ml
Economic model	Low Price

[Click Below](#)

For General Use

Manual pipettes
Pip. Filler/ Stepper
Electronic Dispenser
Magnetic Stirrers
Shakers

For Analytical Labs

Electronic Titrator
Manual Dispenser
Rotary Evaporator
Hotplate
Overhead Stirrers

For Life Science Labs

Electronic Pipettes
SafeVac
Centrifuges-LS
Mixers / Vortex
Thermo Mixers

For Clinical Labs

Centrifuges-Clinical
Rollers and Rotators
Thermo Cyclers
Spectrophotometers



Pipette Filler Levo Plus/Levo ME/Levo

DVAB
Your Accessible Solution

0.45ul Universal Hydrophobic filter



Silicone adapter for 0.1-100ml
pipette tip, fully Autoclavable

Levo ME

Advantages for the user	Features
Easy to operate	Lightweight and single-handed operation
High applicability	For most plastic and glass pipettes from 0.1 to 100ml
Rechargeable battery	Powerful rechargeable Li-ion battery enables 8h runtime

Ordering Information

Levo Plus

Cat. No.	Description
7033100100	Levo Plus with AC adapter, Euro plug
	Levo Plus with AC adapter, North-America plug
	Levo Plus with AC adapter, UK plug

Levo Plus Accessories

Cat. No.	Description
17400032	Levo Plus charger, Euro plug
17000205	Levo Plus Lithium battery
17000103	Filter 0.45um , 1pc/pk (for both Levoplus& Levo ME)

Levo ME

Cat. No.	Description
7013100200	Levo ME with AC adapter, Euro plug
	Levo ME with AC adapter, North-America plug
	Levo ME with AC adapter, UK plug

Levo

Cat. No.	Description
7033201000	Levo - Blue
7033202000	Levo - Red
7033203000	Levo - Green
7033204000	Levo - Yellow

Levo Accessories

Cat. No.	Description
17000110	Filter 3.0µm, 1pc/pk