



## MS-PA

### LED Digital Magnetic Stirrer

#### Features:

- Digital speed control with a range of 100~1500rpm
- Max. stirring quantity of H<sub>2</sub>O at 3L
- LED display shows speeds
- ABS housing provides a good performance of weak acids and alkalis resistance

## MS-PB

### Magnetic Stirrer

#### Features:

- A wide range of 0~1500rpm
- Max. stirring quantity of H<sub>2</sub>O at 3L
- ABS housing provides a good performance of weak acids and alkalis resistance



Chemical resistance



Strong magnetism



LED display



#### Specifications

	MS-PA	MS-PB
Work plate Dimension	φ135(5 inch)mm	φ135(5 inch)mm
Work plate material	ABS	ABS
Motor type	DC motor	DC motor
Motor input power	5W	5W
Motor output power	3W	3W
Power	15W	10W
Voltage	100~120/200~240V,50/60Hz	100~120/200~240V,50/60Hz
Stirring positions	1	1
Max. stirring quantity[H <sub>2</sub> O]	3L	3L
Max. magnetic bar[length]	50mm	50mm
Speed range	100~1500rpm	0~1500rpm
Speed display	LED	Scale
Protection class	IP42	IP42
Dimension [W x D x H]	150×260×80mm	150×260×80 mm
Weight	1.8kg	1.8kg
Permissible ambient temperature and humidity	5~40°C,80%	5~40°C,80%