



Powerful basic model
Delivery pump for small delivery volumes
and high viscosities

Programmable on pump module
 Flexible, thanks to replaceable single and multichannel pump heads

rotarus® standard

| rotarus® standard 50 | Order No. 950 13 22 |
|---|---|
| Characteristics | |
| Flow rate range, depending on the pump head used* | 0,0001 - 1500 ml/min |
| Changeable pump heads | yes |
| Channels | 1 - 12 (Depending on the pump head type used) |
| Range of speed | 0.2 - 100 1/min |
| Operation | TFT touchscreen |
| Program memory | No |
| Pump direction change | yes |
| Speed control | Load-independent |
| Speed configuration | TFT touchscreen |
| Speed configuration (increments)1/min | 0.1 |
| Operation modes | |
| | Pumps according to speed |
| Special function | |
| | Stackable |
| | Usable on stand |
| | Configurable via menu |
| Interfaces | |
| | Analog input for speed specification |
| | Analog output for speed |
| | Digital input start / stop ; rotation direction |
| | Selectable digital output |
| General data | |
| Housing | White powder coated |
| Continuous mode | Yes |
| Motor capacity | 50 Watt |
| Motor type | Stepper motor |
| Max. Power consumption | 60 VA |
| Power supply | 100-240 VAC ; 50-60 Hz |
| Protection class | IP54 |
| Standards | CE |
| Depth/width/height in mm | 205/100/220 |
| Weight | 7.1 kg |

*Indicated flow rates are approximate values measured with water, at 0 meters suction lift and free outlet.