



Order number 51343153

former order number -

Specifications

pH range 1 - 11

zero point (pH) 7.00 ± 0.25

temperature range 0 - 80 °C

temperature probe -

shaft material glass

membrane glass LoT

membrane resistance < 250 MΩ

reference system ARGENTHAL™

reference diaphragm open junction

reference electrolyte XEROLYT®PLUS polymer

bridge electrolyte -

electrode head S7

cable length -

connector -

cell constant (cm⁻¹) -

cell type -

indicating metal -

Storage FRISCOLYT-B®
(51340053)

Special features

high performance puncture electrode:

- XEROLYT®PLUS polymer electrolyte
- open reference junction
- ARGENTHAL™ reference system

Choice of cables

This electrode has a grey S7 (InLab) connector. The cable selector guide offers a large choice of grey S7 cables.

The grey S7 connectors are compatible with the red S7 connectors. The color is only design-related.

replaces the

InLab®427