



**Order number** 51340335  
former order number -

**Specifications**  
pH range -  
zero point (pH) -

temperature range 0 - 80 °C  
temperature probe NTC 30 kΩ

shaft material polycarbonate  
membrane glass -  
membrane resistance -

reference system -  
reference diaphragm -  
reference electrolyte -

bridge electrolyte -

electrode head -  
cable length 1.0 m  
connector Mini-DIN

cell constant (cm<sup>-1</sup>) 0.55  
cell type 4 cell graphite

indicating metal -

**Storage** dry

**Additional specifications**  
measuring range 0.01 - 200 mS/cm