




The **installation drawing** is available under [www.kirsch-medical.com](http://www.kirsch-medical.com) for each model in the download area.



#### **Suitable for integration**

Compact freezer designed for the integration into furniture. Can also be used as a standalone device.



#### **Variable stock options**

due to standard drawer with length and cross dividers and shelves.

# LABO-85



- **External housing** made from galvanised sheet steel (rust-proof), with white, anti-scratch powder coating. Length of the plug cable: approx. 2.8 m.
- **Interior** made from impact-resistant, white plastic with moulded-in shelves.
- **Interior consists** of one drawer with stops and two plastic-coated shelves. There is one adjustable length divider and six adjustable cross dividers per drawer. Storage surface per drawer = 0.11 m<sup>2</sup>.
- **Extra-thick 42 mm energy saving insulation**, made from high-quality, compression-moulded and environmentally friendly material.
- **Lockable door** with easy-to-replace plastic magnetic seal frame.

- **Door hinge** on the right-hand side (by default, see illustration), or the left-hand side, can be retrofitted.
- **Automatic defrosting** with time and temperature monitoring.
- **Condensate evaporates** in the refrigerating machine interior.
- **Mechanical controller** enables variable temperature control. The temperature is maintained automatically, irrespective of changing external temperatures, provided that the ambient temperature is at least approx. 3 °C above the set interior temperature.
- **Statically ventilated refrigerating machine**, hermetically sealed, energy saving, low-noise, easy to service.

## Specification

Capacity	80 litres
Temperature setting	approx. +2 °C to +12 °C
Voltage	220 - 240 V, 50 Hz
Power consumption	100 watts
Normal consumption	0.60 kWh/24 h
Admissible ambient temperature	from +10 °C to +38 °C
Heat emission (max.)	107 watts
Overall dimensions (including distance from wall)	w x d x h = 47 x 51 x 74 cm
Interior dimensions	w x d x h = 38 x 41 x 63 cm (Usable depth above: 5 cm, below: 15 cm less)
Exterior dimensions with door open at 90°	w x d = 47 x 96 cm
Shelf size	w x d = 23 x 37 cm
Drawer inner dimensions	w x d x h = 34 x 32 x 5.6 cm
Max. load drawer/shelf	9.5 kg/25 kg
Weight	Net weight 31 kg, gross weight 37 kg

## Optional equipment

- **Decorative door frame** for attaching the customer's decorative plates (exterior dimensions: 48cm in width).
- **Door fitting.**
- **Equipped with castors.**
- **Additional plastic drawer.**
- **Additional length and cross dividers.**
- **Top drawer can be locked.**
- **Condensate container** can be emptied manually if no condensation is required (e.g., in an operating theatre).
- **Refrigerating machine fan** if static ventilation is not sufficient after installation, or if 60 Hz operation is intended.
- **GSM-Module** for alarm text message transmission (e.g., to a mobile phone).
- **Temperature documentation:** see detailed information on page 56.

